

## Supplemental Material

### Exposome Characterization Suspect-screening Analysis of of Neurotoxic Substance Environmental Chemicals s in Paired Human Cerebrospinal Fluid and Serum Samples

Long Zhang, Pan Yang, Yaqing Shu, Wei Huang, Wenwen Sun, Xiaotu Liu, Da Chen

Formatted: Font: Not Bold

Formatted: Centered

#### Chemicals and reagents

Reference standards of the 28 target environmental chemicals (Table S1), compound standards including 13 per- and polyfluoroalkyl substances (i.e., PFBA, PFHxA, PFHpA, PFOA, PFNA, PFDA, PFUdA, PFDeA, PFBS, PFHxS, PFHpS, PFOS, 6:2 Cl-PFESA), three organophosphate esters (i.e., TBP, TPHP, TCIPP), four personal care products (i.e., MeP, EtP, BuP, TCS), two photoinitiators (i.e., MK, PI-654), three bisphenols (i.e., BPA, BPF, BPS), two phthalate ester metabolites (i.e., MMP, M(i)BP) and one antioxidant (i.e., BHT). The chemicals mentioned above were purchased from Wellington Laboratories (Guelph, Ontario, Canada), AccuStandard (New Haven, Connecticut, USA), and Toronto Research Chemicals (Toronto, Canada), respectively. Twenty-one isotopically labeled reference standards were used as internal standards and purchased from Wellington Laboratories (Guelph, Ontario, Canada), Dr. Ehrenstorfer, LGC Standards (Middlesex, UK), and TRC: Toronto Research Chemicals (Toronto, Canada) (Table S1). Methanol (HPLC grade), ethyl acetate (HPLC grade), hexane (HPLC grade), and water were all high-performance liquid chromatography grade and (UPLC grade) were purchased from Fisher Scientific (Pittsburgh, PA, USA). The sheep serum used for quality control was purchased from Future Biotechnology Company (Guangzhou, China).

Formatted: Not Highlight

Formatted: Not Highlight

#### Sample collection

Our study recruited 180 outpatients diagnosed with mental illness (n = 43), lumbar disc herniation (n = 48), spinal stenosis (n = 37), or viral meningitis (n = 52) from the Third Affiliated Hospital of Sun Yat-Sen University (Guangzhou, China) in 2022–2023. Basic demographic data of the participants are summarized in Table S2. The cerebrospinal fluid (This study was conducted on 180 volunteers recruited at the Department of Neurosurgery, the Third

Affiliated Hospital of Sun Yat-sen University, Guangzhou, China, from April 2022 to May 2023. CSF samples were collected by a senior physician using a clinical lumbar puncture method. Paired blood was also collected from the participants by a specialized nurse in the same day as the CSF collection, and paired serum samples were collected by a specialized nurse. Serum was collected after centrifugation of the blood samples. Field blanks were prepared by collecting deionized water with the same procedures. Collected A total of 180 serum-cerebrospinal fluid pairs were collected, and all samples were stored at  $-80^{\circ}\text{C}$  for subsequent analysis. These samples were obtained from patients with different diagnostic symptoms (mental illness ( $n = 43$ ), lumbar disc herniation ( $n = 48$ ), spinal stenosis ( $n = 37$ ), or viral meningitis ( $n = 52$ )). The ages of the participants were  $39 \pm 19$  years old. All participants completed an informed consent form. The study protocol was approved by the Clinical Research Ethics Committee of the Third Affiliated Hospital of Sun Yat-sen University.

#### Sample pretreatment

An aliquot of 0.2 mL of serum/CSF was spiked with 213 internal standards (with a concentration ranging from 0.05 to 1.4 ng/mL, 5 ng each) and extracted with was performed using 3 mL of a mixture of ethyl acetate and n-hexane (3:2, v/v) containing 0.6% formic acid. The extraction was conducted under sonication for 5 min. After the supernatant was collected, the extraction was repeated twice with the same approach. The collected supernatants after three cycles of extraction were combined. Prior to extraction, isotopically labeled replacement standards (5 ng each) were added to 200  $\mu\text{L}$  of serum or CSF and used to correct for analyte loss during extraction. After extraction, the combined extracts were concentrated to near dryness, and reconstituted with 50  $\mu\text{L}$  of methanol. The extract was frozen overnight at  $-80^{\circ}\text{C}$ , followed by and then centrifugation at  $(15000 \text{ rpm}, -10^{\circ}\text{C})$  for 5 min. The supernatant was retained and reconstituted and an internal standard (5 ng each) was added to destabilize instrument performance by eliminating the effect of inter-sample variability in the sample volume on the instrument response and monitoring. The samples were then diluted to 50  $\mu\text{L}$  and for instrumental analysis was performed.

#### Compound identification with high-resolution mass spectrometer

A single pooled extract was prepared by combining 5  $\mu\text{L}$  each of the 180 CSF or serum extracts. Compound identification in pooled extracts was conducted with a Vanquish ultra-high performance liquid chromatography (UHPLC) interfaced with a ThermoFisher Orbitrap 240 with a H-ESI ion source (Thermo Scientific, Waltham, MA, USA). The UHPLC was equipped

with a Waters ACQUITY UPLC BEH shield RP 18 column (2.1 mm × 100 mm, 1.7 μm particle size; Waters, Milford, USA). The column oven temperature and flow rate were set at 40 °C and 0.4 mL/min, respectively. The mobile phase for positive mode analysis consisted of water containing 0.1% formic acid (A) and methanol (B). For negative mode, the mobile phases consisted of water containing 0.2 mM ammonium acetate (A) and methanol (B). Both modes shared the same gradient, and the elution procedure in positive or negative modes was set as follows: 0–2 min 30% B; 2–7 min, 30% B ramped to 60% B; 7–23 min, 60% B ramped to 100% B; 23–27 min, 100% B; 27–27.1 min, 100% B decreased to 30% B; and 27.1–30 min, 30% B. 30% B (held for 2 min) and linearly ramped to 60% B in 5 min, followed by a linear increase to 100% B in 16 min and held for 4 min, and then changed to 30% B in 0.1 min and equilibrated for 2.9 min. The total run time was 30 min. The column temperature was maintained at 40 °C, LC flow rate was 0.4 mL/min. Mass spectra analysis was conducted in data-dependent MS/MS spectra acquisition method. The ion source conditions were set as following: ion source type, H-ESI; spray voltage, 3500 v for positive mode, 28500 v for negative mode; sheath gas flow rate, 450 arbitrary units; aux gas flow rate, 10 arbitrary units; sweep gas flow rate, 1 arbitrary unit; ion transfer tube temp, 3205°C; vaporizer temperature, 350°C. The following acquisition parameters were used for MS1 analysis: resolution, 60,000, 120,000; scan range 10050–1,700,500 m/z; RF lens (%), 70; AGC target, standard; maximum injection time mode, auto; spectrum data type, profile. Data-dependent MS/MS parameters included: collision energy type, normalized; collision energy mode, stepped; HCD collision energies (%), 20, 40, 100, 60; resolution, 45,000, 30,000; maximum injection time mode, auto; AGC target, standard; number of dependents scans, 10; isolation window (m/z), 1.52; spectrum data type, profile; dynamic exclusion, auto, custom. Ion source type was H-ESI. In DDA mode, the spray voltage was set to 3500 V and 2800 V for positive and negative ionization modes, respectively. The full scan range was 100–1000 Da, the resolution was 120,000. Add database to MassList, filter type was TargetedInclusionMassList, Mass tolerance <5 ppm. Use dynamic exclusions, Mass tolerance <5 ppm. MS/MS Collision Energy Type was Normalized, the resolution was 30,000, fragment ion spectra (MS/MS) were collected at 10, 40, and 100 eV collision energies.

### Quantitative analysis of target compounds

Determination of 28 target compounds was asere conducted on a Shimadzu UHPLC coupled to an AB Sciex 7500 Q Trap triple quadrupole mass spectrometer MS/MS (AB Sciex, Toronto, Canada). Chromatographic separation was achieved by using a Luna Omega 3 μm PS C18 (2) 100 Å column (100 mm × 3.0 mm; 3 μm particle size; Phenomenex, Torrance, CA). The mobile

Formatted: Font: (Default) Times New Roman

Formatted: Font: (Default) Times New Roman

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Font: (Default) Times New Roman

Formatted: Font: (Default) Times New Roman, Not Highlight

Formatted: Not Highlight

Formatted: Font: (Default) Times New Roman

phase consisted of ~~phase A of~~ 0.2 mM ammonium acetate (A) and B-methanol (B). The flow rate was 400 µL/min, and the column temperature was kept at 40 °C. The following gradient was employed: 30% B (held for 2 min) and linearly ramped to 60% B in 5 min, followed by a linear increase to 100% B in 16 min and held for 4 min, and then changed to 30% B in 0.1 min and equilibrated for 2.9 min. The total run time was 30 min. ~~The MS was equipped with a TurboIonSpray® electrospray ionization (ESI) probe and operated in multiple reaction monitoring (MRM) mode (positive and negative). The MS was operated in multiple reaction monitoring (MRM) mode (positive and negative).~~ The ESI source conditions were as follows: ion source Gas 1, 35 psi; ion source Gas 2, 70 psi; curtain gas, 40 psi; collision gas, 9 psi; source temperature, 320 °C; Ion-Spray voltage floating, 2000 V.

Formatted: Font: (Default) Times New Roman

#### Quality assurance and control (QA/QC) procedures

Multiple procedures were performed to ensure data quality, including ~~recovery test, matrix spiking analysis, matrix effect evaluation, monitoring of~~ procedural blank ~~contamination, and assessment of inter-batch variation~~ assessment, and substitution standard recovery monitoring. ~~Considering the rarity of CSF samples, QA/QC tests were conducted only with sheep serum (Guangzhou Future Biotechnology Company, China) as the substitutive matrix for human serum and CSF.~~

Spiking tests were conducted to evaluate the recovery efficiency from sample treatments. An aliquot of 0.2 mL of sheep serum was spiked with 28 target analytes (5 ng each) and internal standards (5 ng each) and processed in five replicates with the aforementioned methodology, along with two matrix blanks (sheep serum spiked with internal standards only). The recovery efficiency (RE) was determined as following.

$$RE (\%) = 100 \times \frac{C_{matrix} - C_{mbk}}{C_{std}} \quad (\text{Eq. 1})$$

Where  $C_{matrix}$ ,  $C_{mbk}$  and  $C_{std}$  represent the determined abundances of target chemicals in sheep plasma (spiked with target analytes), matrix blanks, and the solution of standard mixture, representatively. The recoveries of individual compounds were determined to range from 72 (± 15.5%) to 128.19 (± 5.6%).

Matrix effects were evaluated for the target analytes in serum. In brief, sheep serum (with no standards spiked) was extracted with the aforementioned approach, and each extract was divided into two sub-samples with equal volume. Sub-sample A was spiked with 100 µL of a standard mixture of analytes and sub-sample B was spiked with 100 µL of methanol. An

Formatted: Font: 小四

Formatted: Font: 小四

Formatted: Line spacing: 1.5 lines

Formatted: Font: (Default) Times New Roman

external standard solution (S) was prepared by mixing the 100 µL of analyte mixtures with 100 µL methanol. By comparing the response differences of the analytes in the sub-samples A and B to the responses of the analytes in the external standard, a matrix effect (ME) value was calculated as:

$$ME (\%) = 100 \times \frac{A_i - B_i}{S_i} \quad (\text{Eq. 2})$$

Where  $A_i$ ,  $B_i$  and  $S_i$  represent the chromatographic peak areas of the analyte (i) in sub-samples A and B and external standard solution (S), respectively. Matrix effect was determined to be 74 (± 5XX)% to 119 (± 7.9XX)% (Table SX3).

A laboratory procedural blank was processed along with every 10 samples. A field blank was prepared for each batch of 50 samples by collecting high-performance liquid chromatography grade water with the similar procedures as that for serum/CSF collection. Only TBP, TPHP, TCIPP, MeP, M(i)BP, and BHT were detectable in the procedural or field blanks, with an average concentration ranging from 0.01 (±0.009) to 0.13 7(±0.17) ng/mL. Concentrations of these compounds in serum/CSF samples were reported after subtracting the average contamination levels in the blanks. A QC sample (sheep serum spiked with known amounts of target compounds along with internal standards) was processed along with each batch of 50 samples. The inter-batch coefficients of variation in the recovery of target analytes ranged from -xx2.0% to xx17.0%.

Formatted: Font: (Default) Times New Roman

Formatted: Font: 小四

Formatted: Font: 小四

Formatted: Font: (Default) Times New Roman, (Asian) +Body Asian (等线)

Formatted: Font: (Asian) +Body Asian (等线)

Formatted: Font: (Default) Times New Roman, (Asian) +Body Asian (等线)

Formatted: Indent: First line: 0.74 cm, Don't adjust right indent when grid is defined, Line spacing: 1.5 lines, Don't snap to grid

Formatted: Not Highlight

Formatted: Not Highlight

Formatted: Not Highlight

Formatted: Font: (Default) Times New Roman

Formatted: Line spacing: 1.5 lines

Table S1. Target chemical and corresponding internal standard, including its full name, abbreviation or commercial name, chemical formula, CAS number, and reference standard source. Five laboratory procedure blanks and one matrix spiked sample (sheep serum) were processed for each batch of 90 samples. Only TBP, TPHP, TCIPP, MeP, M(i)BP and BHT were detected in the procedure blanks, with an average concentration range of 0.01–0.17 ng/ml. For these chemicals, their concentrations reported in the article were corrected for blanks. The relative standard deviation of the measured values for individual compounds added to sheep serum was less than 20%. The range of recoveries for the replacement standards in the sample analysis was 72%–128%. The limit of detection (LOD) of an analyte was defined as its response three times the standard deviation of the noise. The limit of detection (LOD) for analytes not present in the program blank was defined as three times the standard deviation of the noise of the response value; the LODs determined ranged from 0.0003–0.34 ng/mL, respectively.

Chemical	Full name	Chemical formula	CAS number	Supplier <sup>a</sup>
<i>Per- and polyfluoroalkyl substances</i>				
PFBA	perfluorobutanoic acid	C <sub>4</sub> HF <sub>7</sub> O <sub>2</sub>	375-22-4	Wellington
PFHxA	perfluorohexanoic acid	C <sub>6</sub> HF <sub>11</sub> O <sub>2</sub>	307-24-4	Wellington
PFHpA	perfluoroheptanoic acid	C <sub>7</sub> HF <sub>13</sub> O <sub>2</sub>	375-85-9	Wellington
PFOA	perfluorooctanoic acid	C <sub>8</sub> HF <sub>15</sub> O <sub>2</sub>	335-67-1	Wellington
PFNA	perfluorononanoic acid	C <sub>9</sub> HF <sub>17</sub> O <sub>2</sub>	375-95-1	Wellington
PFDA	perfluorodecanoic acid	C <sub>10</sub> HF <sub>19</sub> O <sub>2</sub>	335-76-2	Wellington
PFUdA	perfluoroundecanoic acid	C <sub>11</sub> HF <sub>21</sub> O <sub>2</sub>	2058-94-8	Wellington
PFDoA	perfluorododecanoic acid	C <sub>12</sub> HF <sub>23</sub> O <sub>2</sub>	307-55-1	Wellington
PFBS	perfluorobutanesulfonic acid	C <sub>4</sub> HF <sub>9</sub> O <sub>3</sub> S	375-73-5	Wellington
PFHxS	perfluorohexanesulfonic acid	C <sub>6</sub> HF <sub>13</sub> O <sub>3</sub> S	355-46-4	Wellington
PFHpS	perfluoroheptanesulfonic acid	C <sub>7</sub> HF <sub>15</sub> O <sub>3</sub> S	375-92-8	Wellington
PFOS	perfluorooctanesulfonic acid	C <sub>8</sub> HF <sub>17</sub> O <sub>3</sub> S	1763-23-1	Wellington
6:2 Cl-PFESA	6:2 chlorinated perfluoroalkyl ether sulfonic acid	C <sub>8</sub> F <sub>16</sub> SO <sub>4</sub> HCl	756426-58-1	Wellington
<i>Organophosphate esters</i>				
TBP	tributyl phosphate	C <sub>12</sub> H <sub>27</sub> O <sub>4</sub> P	126-73-8	AccuStandard
TPHP	triphenyl phosphate	C <sub>18</sub> H <sub>15</sub> O <sub>4</sub> P	115-86-6	AccuStandard
TCIPP	tris(2-chloroisopropyl) phosphate	C <sub>6</sub> H <sub>18</sub> Cl <sub>3</sub> O <sub>4</sub> P	13674-84-5	AccuStandard
<i>Personal care products</i>				
MeP	methyl paraben	C <sub>8</sub> H <sub>8</sub> O <sub>3</sub>	99-76-3	AccuStandard
EtP	ethyl paraben	C <sub>9</sub> H <sub>10</sub> O <sub>3</sub>	120-47-8	AccuStandard
BuP	butyl paraben	C <sub>11</sub> H <sub>14</sub> O <sub>3</sub>	94-26-8	AccuStandard
TCS	triclosan	C <sub>12</sub> H <sub>7</sub> Cl <sub>5</sub> O <sub>2</sub>	3380-34-5	AccuStandard

Photoinitiators

MK	4,4'-bis(dimethylamino)benzophenone	C <sub>17</sub> H <sub>20</sub> N <sub>2</sub> O	90-94-8	TRC
PI-651	2,2-dimethoxyphenyl acetophenone	C <sub>16</sub> H <sub>16</sub> O <sub>3</sub>	24650-42-8	TRC

Bisphenols

BPA	2,2-bis(4-hydroxyphenyl)propane	C <sub>15</sub> H <sub>16</sub> O <sub>2</sub>	80-05-7	AccuStandard
BPF	4,4'-dihydroxydiphenylmethane	C <sub>13</sub> H <sub>12</sub> O <sub>2</sub>	620-92-8	AccuStandard
BPS	4,4'-dihydroxydiphenylsulfone	C <sub>12</sub> H <sub>10</sub> O <sub>2</sub> S	80-09-1	AccuStandard

Phthalate ester metabolites

MMP	monomethyl phthalate	C <sub>8</sub> H <sub>8</sub> O <sub>4</sub>	4376-18-5	AccuStandard
M(i)BP	mono-(iso)butyl phthalate	C <sub>12</sub> H <sub>14</sub> O <sub>4</sub>	30833-53-5	AccuStandard

Antioxidants

BHT	2,6-di-tert-butyl-4-hydroxytoluene	C <sub>15</sub> H <sub>24</sub> O	128-37-0	AccuStandard
-----	------------------------------------	-----------------------------------	----------	--------------

Internal standards

MPFBA	perfluoro[(13)C4]butanoic acid	[ <sup>13</sup> C <sub>4</sub> ]HF <sub>7</sub> O <sub>2</sub>	NA <sup>b</sup>	Wellington
MPFHxA	perfluoro[1,2-(13)C2]hexanoic acid	[ <sup>13</sup> C <sub>2</sub> ]C <sub>6</sub> HF <sub>11</sub> O <sub>2</sub>	NA	Wellington
MPFHxS	perfluorohexane[(18)O2]sulfonic acid	C <sub>6</sub> HF <sub>13</sub> [ <sup>18</sup> O <sub>2</sub> ]OS	1585941-14-5	Wellington
MPFOA	perfluoro[1,2,3,4-(13)C4]octanoic acid	[ <sup>13</sup> C <sub>4</sub> ]C <sub>8</sub> HF <sub>15</sub> O <sub>2</sub>	960315-48-4	Wellington
MPFNA	perfluoro[1,2,3,4,5-(13)C5]nonanoic acid	[ <sup>13</sup> C <sub>5</sub> ]C <sub>9</sub> HF <sub>17</sub> O <sub>2</sub>	NA	Wellington
MPFDA	perfluoro[1,2-(13)C2]decanoic acid	[ <sup>13</sup> C <sub>2</sub> ]C <sub>10</sub> HF <sub>19</sub> O <sub>2</sub>	NA	Wellington
MPFUdA	perfluoro[1,2-(13)C2]undecanoic acid	[ <sup>13</sup> C <sub>2</sub> ]C <sub>11</sub> HF <sub>21</sub> O <sub>2</sub>	960315-51-9	Wellington
MPFDoA	perfluoro[1,2-(13)C2]dodecanoic acid	[ <sup>13</sup> C <sub>2</sub> ]C <sub>12</sub> HF <sub>23</sub> O <sub>2</sub>	NA	Wellington
MPFOS	perfluoro[1,2,3,4-(13)C4]octanesulfonic acid	[ <sup>13</sup> C <sub>4</sub> ]C <sub>8</sub> HF <sub>17</sub> O <sub>3</sub> S	960315-53-1	Wellington
TBP-d <sub>27</sub>	tri-n-butyl phosphate-d <sub>27</sub>	C <sub>12</sub> [ <sup>2</sup> H <sub>27</sub> ]O <sub>4</sub> P	61196-26-7	Wellington
TPHP-d <sub>15</sub>	triphenyl phosphate-d <sub>15</sub>	C <sub>18</sub> [ <sup>2</sup> H <sub>15</sub> ]O <sub>4</sub> P	1173020-30-8	Wellington
TDCIPP-d <sub>15</sub>	tris(1,3-dichloro-2-propyl) phosphate-d <sub>15</sub>	C <sub>9</sub> [ <sup>2</sup> H <sub>15</sub> ]Cl <sub>6</sub> O <sub>4</sub> P	1447569-77-8	Wellington
BHT-d <sub>21</sub>	2,6-di(tert-butyl)-4-methyl-phenol-d <sub>21</sub>	C <sub>15</sub> H <sub>23</sub> [ <sup>2</sup> H <sub>21</sub> ]O	64502-99-4	Dr. Ehrenstorfer
BPA-d <sub>6</sub>	bisphenol A-d <sub>6</sub>	C <sub>15</sub> H <sub>10</sub> [ <sup>2</sup> H <sub>6</sub> ]O <sub>2</sub>	86588-58-1	TRC
BPF-d <sub>10</sub>	bisphenol F-d <sub>10</sub>	C <sub>12</sub> H <sub>8</sub> [ <sup>2</sup> H <sub>10</sub> ]O <sub>2</sub>	1794786-93-8	TRC
BPS-13C <sub>12</sub>	bisphenol S- <sup>13</sup> C <sub>12</sub>	[ <sup>13</sup> C <sub>12</sub> ]H10O <sub>2</sub> S	1991267-29-8	TRC
MEPA-d <sub>4</sub>	monoethyl phthalate-d <sub>4</sub>	C <sub>10</sub> H <sub>8</sub> [ <sup>2</sup> H <sub>4</sub> ]O <sub>4</sub>	1219806-03-7	TRC
MeP-d <sub>4</sub>	methyl paraben-d <sub>4</sub>	C <sub>8</sub> H <sub>8</sub> [ <sup>2</sup> H <sub>4</sub> ]O <sub>3</sub>	362049-51-2	TRC
EtP-d <sub>4</sub>	ethyl paraben-d <sub>4</sub>	C <sub>9</sub> H <sub>10</sub> [ <sup>2</sup> H <sub>4</sub> ]O <sub>3</sub>	1219795-53-5	TRC
Bup-d <sub>6</sub>	butyl-paraben-d <sub>6</sub>	C <sub>11</sub> H <sub>12</sub> [ <sup>2</sup> H <sub>6</sub> ]O <sub>3</sub>	1216904-65-2	TRC
TCS-d <sub>3</sub>	triclosan-d <sub>3</sub>	C <sub>12</sub> H <sub>7</sub> [ <sup>2</sup> H <sub>3</sub> ]Cl <sub>3</sub> O <sub>2</sub>	1020719-98-5	TRC

Internal standards

<u>M8PFOA</u>	<u>perfluoro[<sup>13</sup>C<sub>8</sub>]octanoic acid</u>	<u>[<sup>13</sup>C<sub>8</sub>]HF<sub>15</sub>O<sub>3</sub></u>	<u>NA</u>	<u>Wellington</u>
<u>BPA-d<sub>16</sub></u>	<u>bisphenol A-d<sub>16</sub></u>	<u>C<sub>15</sub>[<sup>2</sup>H<sub>16</sub>]O<sub>3</sub></u>	<u>96210-87-6</u>	<u>Sigma-Aldrich</u>
<u>CMP-d<sub>10</sub></u>	<u>coumaphos-d<sub>10</sub></u>	<u>C<sub>14</sub>H<sub>6</sub>[<sup>2</sup>H<sub>10</sub>]ClO<sub>6</sub>PS</u>	<u>287397-86-8</u>	<u>TRC</u>

<sup>3</sup>Suppliers: Wellington: Wellington Laboratories (Guelph, Ontario, Canada); AccuStandard: AccuStandard Inc. (New Haven, Connecticut, USA); TRC: [TRC](#); Toronto Research Chemicals (Toronto, Canada); Dr. Ehrenstorfer: LGC Standards (Middlesex, UK); Sigma-Aldrich: Sigma-Aldrich, Inc. (St. Louis, Missouri, USA). <sup>4</sup>NA: not available.

- Formatted: Font: 10 pt
- Formatted
- Formatted
- Formatted: Line spacing: single
- Formatted: Font: 10 pt
- Formatted: Font: 10 pt, Subscript
- Formatted
- Formatted
- Formatted: Line spacing: single
- Formatted: Font: 10 pt
- Formatted: Font: 10 pt, Subscript
- Formatted
- Formatted
- Formatted: Line spacing: single
- Formatted: Font: 10 pt
- Formatted: Line spacing: single, Don't adjust space between Latin and Asian text, Don't adjust space between Asian text and numbers, Snap to grid
- Formatted: Font: 11 pt, Italic, Superscript
- Formatted: Font: 11 pt
- Formatted: Line spacing: Multiple 1.15 li
- Formatted



Table S2. Demographic characteristics of study participants.<sup>a</sup>

variables	n (%)	Mean ± SD	range
age (years)	180	39 ± 19	13-85
gender			
female	95 (52.8)		
male	85 (47.2)		
disease type			
mental illness	43 (23.9)		
lumbar disc herniation	48 (26.7)		
spinal stenosis	37 (20.6)		
or viral meningitis	52 (28.9)		

<sup>a</sup>Continuous variables are presented as mean ± SD. Category variables are presented as numbers (percentage).

Formatted: Font: 小四

Formatted: Font: Italic, Superscript

Formatted: Font: 小四

Formatted: Font: 10 pt

Formatted: Line spacing: Multiple 1.15 li

Formatted: Centered, Line spacing: single

Formatted Table

Formatted: Font: 10 pt

Formatted: Font: 10 pt

Formatted: Left, Line spacing: single

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Left, Line spacing: single

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Left, Line spacing: single

Formatted: Centered, Line spacing: single

Formatted Table

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 10 pt

Formatted: Centered, Line spacing: single

Formatted: Font: 11 pt, Italic, Superscript

Formatted: Font: 11 pt

Formatted: Line spacing: Multiple 1.15 li

Table S3. Quality assurance and control test results. Method validation results of the chemicals.

Chemical	Blank <sup>a</sup> (ng/mL, n = 18)	Recovery (%, n = 8)	Matrix effect (%, n=10)	Precision	
				Intra-batch (%, n=4)	Inter-batch (%, n=4)
<i>Per- and polyfluoroalkyl substances</i>					
PFBA		78.0 ± 5.4	100 ± 7.8	6.3	4.0
PFHxA		110 ± 4.4	86.1 ± 3.4	4.0	10.3
PFHpA		75.0 ± 29.1	92.3 ± 6.2	8.9	18.0
PFOA		103 ± 6.4	101 ± 10.3	6.2	11.0
PFNA		96.8 ± 5.3	97.7 ± 8.4	4.5	6.0
PFDA		99.8 ± 7.5	100 ± 3.9	12.1	2.0
PFUdA		108 ± 3.9	83.4 ± 3.8	6.0	2.0
PFDoA		109 ± 10.2	85.5 ± 7.7	5.0	5.2
PFBS		110 ± 10.1	103 ± 5.0	6.2	11.2
PFHxS		72.0 ± 15.0	74.0 ± 5.0	9.0	13.0
PFHpS		104 ± 8.3	101 ± 8.4	4.8	9.2
PFOS		90.5 ± 2.6	81.0 ± 6.2	3.0	7.0
6:2 Cl-PFESA		97.1 ± 8.8	79.5 ± 4.5	6.0	3.0
<i>Organophosphate esters</i>					
TBP	0.01 ± 0.009	128 ± 5.6	107 ± 5.2	5.2	2.0
TPHP	0.03 ± 0.05	80.1 ± 2.6	100 ± 6.7	7.0	2.0
TCIPP	0.04 ± 0.03	74.3 ± 4.0	111 ± 13.7	6.0	2.0
<i>Personal care products</i>					
MeP	0.06 ± 0.05	89.9 ± 12.0	119 ± 7.9	17.2	12.9
EtP		90.9 ± 6.0	108 ± 8.2	10.4	6.0
BuP		79.0 ± 5.9	111 ± 4.4	7.0	9.0
TCS		83.0 ± 9.7	79.4 ± 6.7	15.8	17.0
<i>Photoinitiators</i>					
MK		111 ± 9.4	113 ± 6.8	14.0	12.0
PI-651		102 ± 6.7	91.0 ± 25.7	8.2	5.6
<i>Bisphenols</i>					
BPA		120 ± 16.1	103 ± 14.1	3.0	7.0
BPF		118 ± 14.0	110 ± 5.0	9.1	5.8
BPS		101 ± 3.6	116 ± 8.4	11.0	6.0
<i>Phthalate ester metabolites</i>					
MMP		93.0 ± 8.4	99.3 ± 8.6	6.0	7.0
M(i)BP	0.13 ± 0.17	108 ± 9.9	105 ± 13.0	12.1	10.3
<i>Antioxidants</i>					
BHT	0.11 ± 0.19	82.0 ± 5.8	108 ± 4.6	9.3	15.6

<sup>a</sup>No fill indicates no detectable blank contamination.